

Anti-Clec12A Polyclonal Antibody For Wb, If-P, Elisa - Q5Qgz9

Item Number: CM0029919



Introduction

Rabbit polyclonal antibody against human CLEC12A (CLL-1, MICL), validated for WB, IF-P, ELISA, cross-reacts with mouse and rat. Ideal for myeloid cell activation and immune regulation studies.

[Learn More](#)

Parameter	Value
Product Name	CLEC12A Rabbit pAb
Remarks/Alias	CLL1; MICL; CD371; CLL-1; DCAL-2; CLEC12A
Species	Human
GeneID (Human)	160364
GeneID	160364
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 66-265 of human CLEC12A (NP_612210.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5QGZ9
Protein Weight	31kDa
Shipping	Ice bag