

Anti-Clec1B Polyclonal Antibody For Wb, Elisa - Q9P126

Item Number: CM0029272



Introduction

High-specificity Anti-CLEC1B Polyclonal Antibody validated for WB and ELISA. Detects human, mouse, and rat CLEC1B, a platelet PDPN receptor involved in thrombosis and vascular biology. Ideal for immunodiagnostic and research use.

[Learn More](#)

Parameter	Value
Product Name	CLEC1B Rabbit pAb
Remarks/Alias	CLEC2; CLEC2B; PRO1384; QDED721; 1810061I13Rik; CLEC1B
Species	Human
GeneID	51266
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-196 of human CLEC1B (NP_001092901.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P126
Protein Weight	27kDa
Shipping	Ice bag