


Anti-Clec2D Rabbit Polyclonal Antibody For Wb And Elisa - Q9Uhp7

Item Number: CM0016051



Introduction

Rabbit polyclonal antibody against human CLEC2D (LLT1/OCIL), validated for WB and ELISA. Cross-reacts with mouse. Target is a NK cell receptor involved in immune modulation and osteoclast inhibition. Ideal for immunology and bone research.

[Learn More](#)

Parameter	Value
Product Name	CLEC2D Rabbit pAb
Remarks/Alias	CLAX, LLT1, OCIL, CLEC2D
Species	Human
Gene ID	29121 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-154 of human CLEC2D (NP_037401.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UHP7
Protein Weight	22kDa
Shipping	Ice bag