


Anti-Clec3B/Tetranectin Monoclonal Antibody For Wb, If/Icc, Elisa - P05452

Item Number: CM0008963



Introduction

Rabbit monoclonal antibody against human CLEC3B/Tetranectin (P05452), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studying retinal function, exocytotic packaging, and plasminogen binding.

[Learn More](#)

Parameter	Value
Product Name	CLEC3B/Tetranectin Rabbit mAb
Remarks/Alias	TN; TNA; MCDR4; CLEC3B/Tetranectin
Species	Human
GeneID (Human)	7123
GeneID	7123
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05452
Protein Weight	21kDa
Shipping	Ice bag