


Anti-Clec3B/Tetranectin Polyclonal Antibody For If/Icc, Elisa - P05452

Item Number: CM0025057



Introduction

Rabbit polyclonal antibody against human CLEC3B/Tetranectin, suitable for IF/ICC and ELISA. Targets the secreted C-type lectin involved in plasminogen binding, retinal function, and exocytosis. Ideal for studying plasma protein interactions.

[Learn More](#)

Parameter	Value
Product Name	CLEC3B/Tetranectin Rabbit pAb
Remarks/Alias	TN; TNA; MCDR4; CLEC3B/Tetranectin
Species	Human
GeneID (Human)	7123
GeneID	7123
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-202 of human CLEC3B/Tetranectin (NP_003269.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P05452
Protein Weight	21kDa
Shipping	Ice bag