


Anti-Clec4D Polyclonal Antibody For Wb, If-P, Elisa - Q8Wxi8

Item Number: CM0026887



Introduction

High-quality CLEC4D Rabbit pAb validated for WB, IF-P, ELISA. Recognizes human, mouse, rat CLEC4D (CD368, Dectin-3). Ideal for innate immunity and pattern recognition receptor research.

[Learn More](#)

Parameter	Specification
Product Name	CLEC4D Rabbit pAb
Remarks/Alias	MCL, MPCL, CD368, CLEC6, CLEC-6, CLECSF8, Dectin-3, CLEC4D
Species	Human
GeneID (Human)	338339
GeneID	338339
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 39-215 of human CLEC4D (NP_525126.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WXI8
Protein Weight	25kDa
Shipping	Ice bag