

Anti-Clec4E Polyclonal Antibody For Wb, Elisa - Q9R0Q8

Item Number: CM0027975

Introduction

Rabbit polyclonal antibody targeting mouse CLEC4E (Mincle), a C-type lectin innate immune receptor. Suitable for WB and ELISA, with rat cross-reactivity. Ideal for pattern recognition and host defense studies.

[Learn More](#)

Parameter	Value
Product Name	CLEC4E Rabbit pAb
Remarks/Alias	Mincle; Clecsf9
Species	Mouse
GeneID (Human)	26253
GeneID	56619
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-200 of mouse CLEC4E/MINCLE (NP_064332.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q9R0Q8
Protein Weight	24kDa
Shipping	Ice bag