

Anti-Clec4M Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H2X3

Item Number: CM0014423



Introduction

Rabbit polyclonal antibody targeting human CLEC4M (CD299, L-SIGN, DC-SIGNR). Validated for WB, IF/ICC, ELISA; cross-reacts with Human, Mouse, Rat. Research use for pathogen recognition and immune surveillance studies.

[Learn More](#)

Parameter	Specification
Product Name	CLEC4M Rabbit pAb
Remarks/Alias	CD299; LSIGN; CD209L; L-SIGN; DCSIGNR; HP10347; DC-SIGN2; DC-SIGNR; CLEC4M
Species	Human
Gene ID (Human)	10332
Gene ID	10332
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 71-399 of human CLEC4M (NP_055072.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H2X3
Protein Weight	45kDa
Shipping	Ice bag