


Anti-Clec5A Rabbit Polyclonal Antibody For Wb And Elisa - Q9Ny25

Item Number: CM0028197



Introduction

High-quality rabbit polyclonal antibody against human CLEC5A (MDL-1), validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Ideal for innate immunity, viral infection, and osteoclastogenesis research.

[Learn More](#)

Parameter	Value
Product Name	CLEC5A Rabbit pAb
Remarks/Alias	MDL1; MDL-1; CLECSF5; CLEC5A
Species	Human
Gene ID (Human)	23601
Gene ID	23601
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-188 of human CLEC5A (NP_037384.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NY25
Protein Weight	22kDa
Shipping	Ice bag