

Anti-Clec7A Polyclonal Antibody For Wb, Elisa - Q9Bxn2

Item Number: CM0026675

Introduction

Rabbit polyclonal antibody against human CLEC7A (Dectin-1). Validated for WB and ELISA, cross-reacts with human, mouse, rat. Ideal for studies on innate immunity, fungal recognition, and beta-glucan signaling.

[Learn More](#)

Parameter	Value
Product Name	CLEC7A Rabbit pAb
Remarks/Alias	BGR; CD369; CANDF4; SCARE2; DECTIN1; CLECSF12; CLEC7A
Species	Human
GeneID (Human)	64581
GeneID	64581
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 140-240 of human CLEC7A (NP_922938.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BXN2
Protein Weight	28kDa
Shipping	Ice bag