


Anti-Clc1 Polyclonal Antibody For Wb, If/Icc, Elisa - 000299

Item Number: CM0016992



Introduction

Rabbit polyclonal antibody against human CLIC1 (Chloride intracellular channel protein 1) for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Useful for ion channel and redox biology research.

[Learn More](#)

Parameter	Specification
Product Name	CLIC1 Rabbit pAb
Remarks/Alias	G6; CL1C1; NCC27; CLCNL1; CLIC1
Species	Human
Gene ID (Human)	1192
Gene ID	1192
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00299
Protein Weight	27kDa
Shipping	Ice bag