


Anti-Clc1 Rabbit Polyclonal Antibody For Wb, Elisa - 000299

Item Number: CM0029653



Introduction

Rabbit polyclonal antibody targeting human CLIC1 (27 kDa). Validated for WB and ELISA; cross-reacts with mouse and rat CLIC1. Ideal for ion channel and oxidoreductase research.

[Learn More](#)

Parameter	Value
Product Name	CLIC1 Rabbit pAb
Remarks/Alias	G6; CL1C1; NCC27; CLCNL1; CLIC1
Species	Human
GeneID (Human)	1192
GeneID	1192
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 42-241 of human CLIC1 (NP_001279.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00299
Protein Weight	27kDa
Shipping	Ice bag