

## Anti-Clc4 Rabbit Monoclonal Antibody For Wb And Elisa - Q9Y696

Item Number: CM0002598

### Introduction

Recombinant anti-CLIC4 rabbit monoclonal antibody validated for Western blot (WB) and ELISA. Detects human, mouse, rat CLIC4 (~29 kDa). Ideal for studying ion channel function, redox regulation, and angiogenesis.

[Learn More](#)

| Parameter        | Value                                                                                                                                 |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | CLIC4 Rabbit mAb                                                                                                                      |
| Remarks/Alias    | H1; huH1; p64H1; CLIC4L; MTCLIC; CLIC4                                                                                                |
| Species          | Human                                                                                                                                 |
| Gene ID (Human)  | 25932                                                                                                                                 |
| Gene ID          | 25932                                                                                                                                 |
| Immunogen        | Recombinant protein/Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human CLIC4 (NP_039234.1). |
| Source           | Rabbit                                                                                                                                |
| Category         | Monoclonal Antibodies                                                                                                                 |
| Application      | WB, ELISA                                                                                                                             |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                     |
| SWISS            | Q9Y696                                                                                                                                |
| Protein Weight   | 29kDa                                                                                                                                 |
| Shipping         | Ice bag                                                                                                                               |