


Anti-CLIC4 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y696

Item Number: CM0028391



Introduction

Rabbit polyclonal antibody targeting human CLIC4 (p64H1), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for chloride intracellular channel, oxidoreductase, and angiogenesis research.

[Learn More](#)

Parameter	Value
Product Name	CLIC4 Rabbit pAb
Remarks/Alias	H1; huH1; p64H1; CLIC4L; MTCLIC; CLIC4
Species	Human
GeneID (Human)	25932
GeneID	25932
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-253 of human CLIC4 (NP_039234.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y696
Protein Weight	29kDa
Shipping	Ice bag