


Anti-Clip1 Polyclonal Antibody For Wb, Elisa - P30622

Item Number: CM0017178



Introduction

Human CLIP1 rabbit polyclonal antibody validated for WB and ELISA, cross-reactive with mouse and rat. Targets microtubule plus-end tracking protein involved in cytoskeleton dynamics and vesicle trafficking. Ideal for cancer and cell biology research.

[Learn More](#)

Parameter	Value
Product Name	CLIP1 Rabbit pAb
Remarks/Alias	RSN, CLIP, CYLN1, CLIP170, CLIP-170, CLIP1
Species	Human
GeneID (Human)	6249
GeneID	6249
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1228-1427 of human CLIP1 (NP_002947.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30622
Protein Weight	162kDa
Shipping	Ice bag