


Anti-Clip1 Rabbit Polyclonal Antibody For Wb, Elisa - P30622

Item Number: CM0030379



Introduction

Rabbit polyclonal antibody targeting human CAP-Gly domain-containing linker protein 1 (CLIP1/CLIP-170). Validated for Western blot and ELISA; cross-reacts with mouse and rat. Suitable for microtubule dynamics and vesicle trafficking research.

[Learn More](#)

Parameter	Value
Product Name	CLIP1 Rabbit pAb
Remarks/Alias	RSN, CLIP, CYLN1, CLIP170, CLIP-170, CLIP1
Species	Human
Gene ID (Human)	6249
Gene ID	6249
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1228-1427 of human CLIP1 (NP_002947.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30622
Protein Weight	162kDa
Shipping	Ice bag