

Anti-CLN8 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Uby8

Item Number: CM0021701

Introduction

Rabbit anti-CLN8 polyclonal antibody for WB, IHC-P, ELISA. Targets human CLN8 (Q9UBY8), an ER membrane protein involved in neuroprotection and cell proliferation. Suitable for neuronal ceroid lipofuscinosis research.

[Learn More](#)

Parameter	Value
Product Name	CLN8 Rabbit pAb
Remarks/Alias	EPMR; TLCD6; C8orf61; CLN8
Species	Human
GeneID (Human)	2055
GeneID	2055
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q9UBY8
Protein Weight	33kDa
Shipping	Ice bag