

Anti-Clock Rabbit Monoclonal Antibody For Wb, Elisa - O15516

Item Number: CM0012319



Introduction

Anti-CLOCK rabbit monoclonal antibody validated for WB and ELISA. Reacts with human, mouse, rat CLOCK. Synthetic peptide immunogen (aa 100-200). Ideal for circadian rhythm and transcription regulation research.

[Learn More](#)

Parameter	Value
Product Name	CLOCK Rabbit mAb
Remarks/Alias	KAT13D; bHLHe8; CLOCK
Species	Human
GeneID (Human)	9575
GeneID	9575
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human CLOCK (NP_004889.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15516
Protein Weight	95kDa
Shipping	Ice bag