

Anti-Clock Rabbit Polyclonal Antibody For Wb, Elisa - O15516

Item Number: CM0028150



Introduction

Rabbit polyclonal antibody targeting human CLOCK (KAT13D/bHLHe8), validated for Western blot and ELISA applications. Detects CLOCK in human, mouse, and rat samples. Suitable for circadian rhythm and transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	CLOCK Rabbit pAb
Remarks/Alias	KAT13D; bHLHe8; CLOCK
Species	Human
GeneID (human)	9575
GeneID	9575
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 450-600 of human CLOCK (NP_004889.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15516
Protein Weight	95kDa
Shipping	Ice bag