


Anti-Clpb Rabbit Polyclonal Antibody For Ihc-P, Elisa - Q9H078

Item Number: CM0022714



Introduction

Rabbit polyclonal antibody targeting human CLPB, validated for IHC-P and ELISA. Detects mouse and rat CLPB. Useful for studying mitochondrial protein homeostasis and innate antiviral immunity.

[Learn More](#)

Parameter	Value
Product Name	CLPB Rabbit pAb
Remarks/Alias	CLPB; ANKCLB; HSP78; MEGCANN; MGCA7; SKD3; caseinolytic peptidase B protein homolog
Species	Human
GeneID (Human)	81570
GeneID	81570
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 572-707 of human CLPB (NP_110440.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9H078
Protein Weight	57kDa/72kDa/74kDa/75kDa/78kDa
Shipping	Ice bag