


Anti-Clpp Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q16740

Item Number: CM0009114



Introduction

CLPP Rabbit Monoclonal Antibody for Western blot (WB), IHC-P, IF/ICC, and ELISA. Validated in human, mouse, and rat. Targets mitochondrial ClpXP protease subunit, 30 kDa protein, key for protein quality control.

[Learn More](#)

Parameter	Value
Product Name	CLPP Rabbit mAb
Remarks/Alias	DFNB81; PRLTS3; CLPP
Species	Human
GeneID (Human)	8192
GeneID	8192
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 178-277 of human CLPP (Q16740).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16740
Protein Weight	30kDa
Shipping	Ice bag