


Anti-Clpp Rabbit Polyclonal Antibody For Wb, Elisa - Q16740

Item Number: CM0029763



Introduction

Rabbit polyclonal antibody against human CLPP (mitochondrial ClpP protease) for Western blot and ELISA applications. Cross-reacts with human, mouse, rat. Useful for mitochondrial protein quality control research.

[Learn More](#)

Parameter	Value
Product Name	CLPP Rabbit pAb
Remarks/Alias	DFNB81; PRLTS3; CLPP
Species	Human
GeneID (Human)	8192
GeneID	8192
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16740
Protein Weight	30kDa
Shipping	Ice bag