

Anti-Clps (Colipase) Rabbit Polyclonal Antibody For Wb And Elisa - P04118

Item Number: CM0015995

Introduction

CamelBio provides rabbit polyclonal antibody against CLPS (colipase), validated for WB and ELISA with cross-reactivity to mouse and rat. Suitable for research in lipid metabolism, digestion, and satiety signaling.

[Learn More](#)

Parameter	Value
Product Name	CLPS Rabbit pAb
Remarks/Alias	CLPS; colipase
Species	Human
GenelD (Human)	1208
GenelD	1208
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P04118
Protein Weight	12kDa
Shipping	Ice bag