

Anti-Clptm1 Rabbit Polyclonal Antibody For Wb And Elisa - 096005

Item Number: CM0025129

Introduction

Rabbit polyclonal antibody targeting human CLPTM1. Suitable for Western blot (WB) and ELISA. Cross-reacts with mouse and rat. Use to investigate GABAergic transmission, lipid scrambling, and GPI-anchor processing.

[Learn More](#)

Parameter	Specification
Product Name	CLPTM1 Rabbit pAb
Remarks/Alias	CLPTM1
Species	Human
Gene ID (Human)	1209
Gene ID	1209
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human CLPTM1 (NP_001285.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	096005
Protein Weight	76kDa
Shipping	Ice bag