


Anti-Clpx Rabbit Monoclonal Antibody For Wb, Elisa - 076031

Item Number: CM0003266



Introduction

Rabbit monoclonal antibody against human CLPX, mitochondrial ATP-dependent chaperone. Validated for WB and ELISA, cross-reacts with mouse and rat. Suitable for studying mitochondrial protein quality control and heme biosynthesis.

[Learn More](#)

Parameter	Value
Product Name	CLPX Rabbit mAb
Remarks/Alias	EPP2; CLPX
Species	Human
GeneID (Human)	10845
GeneID	10845
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 534-633 of human CLPX (076031).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	076031
Protein Weight	69 kDa
Shipping	Ice bag