

Anti-Clstn3 Polyclonal Antibody For Wb, Elisa - Q9Bqt9

Item Number: CM0013051

Introduction

High-quality rabbit polyclonal antibody targeting human Calsyntenin-3 (CLSTN3), suitable for Western blot and ELISA applications. Detects mouse CLSTN3, with a predicted molecular weight of 106 kDa. Ideal for synapse biology and thermogenesis studies.

[Learn More](#)

Parameter	Specification
Product Name	CLSTN3 Rabbit pAb
Remarks/Alias	CSTN3; CDHR14; alcbeta; CLSTN3
Species	Human
GeneID (Human)	9746
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 730-850 of human CLSTN3 (NP_055533.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9BQT9
Protein Weight	106kDa
Shipping	Ice bag