

Anti-Clta Polyclonal Antibody For Wb, Elisa - P09496

Item Number: CM0028616



Introduction

Rabbit polyclonal antibody targeting human Clathrin light chain A (CLTA). Validated for WB and ELISA, with cross-reactivity to mouse and rat. Supports endocytosis and mitotic spindle research.

[Learn More](#)

Parameter	Value
Product Name	CLTA Rabbit pAb
Remarks/Alias	LCA; CLTA
Species	Human
GeneID (Human)	1211
GeneID	1211
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-218 of human CLTA (NP_001824.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09496
Protein Weight	27kDa
Shipping	Ice bag