

Anti-Cltb Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P09497

Item Number: CM0029620

Introduction

Rabbit polyclonal antibody against human CLTB (Clathrin light chain B), validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat samples. Suitable for research on coated pits and vesicles.

[Learn More](#)

Parameter	Value
Product Name	CLTB Rabbit pAb
Remarks/Alias	LCB; CLTB
Species	Human
Gene ID (Human)	1212
Gene ID	1212
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-211 of human CLTB (NP_001825.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09497
Protein Weight	25kDa
Shipping	Ice bag