

Anti-Cluh Rabbit Polyclonal Antibody (Wb, Ihc-P, If/Icc, Elisa) - O75153

Item Number: CM0014260

Introduction

Rabbit polyclonal antibody targeting human CLUH (CLU1), suitable for WB, IHC-P, IF/ICC, and ELISA. Detects human and mouse CLUH. Useful for studying mitochondrial biogenesis and cytoplasmic distribution.

[Learn More](#)

Parameter	Value
Product Name	CLUH Rabbit pAb
Remarks/Alias	CLU1; CLUH
Species	Human
GenelD (Human)	23277
GenelD	23277
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1120-1309 of human CLUH (NP_056044.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O75153
Protein Weight	147kDa
Shipping	Ice bag