


Anti-Clusterin Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P10909

Item Number: CM0017974



Introduction

Rabbit polyclonal antibody against Clusterin (ApoJ, TRPM-2), validated for WB, IHC-P, IF/ICC, ELISA; cross-reacts with human, mouse, rat. Ideal for apoptosis, complement, and neurodegeneration studies.

[Learn More](#)

Parameter	Value
Product Name	Clusterin Rabbit pAb
Remarks/Alias	CLI; AAG4; APOJ; CLU1; CLU2; KUB1; SGP2; APO-J; SGP-2; SP-40; TRPM2; TRPM-2; NA1/NA2; Clusterin
Species	Human
GeneID (human)	1191
GeneID	1191
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 220-449 of human Clusterin alpha chain (NP_976084.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10909
Protein Weight	52kDa
Shipping	Ice bag