


Anti-Clusterin Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P10909

Item Number: CM0002171



Introduction

Recombinant rabbit monoclonal anti-Clusterin antibody validated for WB, IHC-P, ELISA. Detects human Clusterin, an extracellular chaperone involved in apoptosis, complement inhibition, and protein aggregation. Suitable for diagnostic and research applications.

[Learn More](#)

Parameter	Details
Product Name	Clusterin Rabbit mAb
Remarks/Alias	CLI, AAG4, APOJ, CLU1, CLU2, KUB1, SGP2, APO-J, SGP-2, SP-40, TRPM2, TRPM-2, NA1/NA2, Clusterin
Species	Human
Gene ID (Human)	1191
Gene ID	1191
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 314-449 of human Clusterin alpha chain (NP_001822.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P10909
Protein Weight	52kDa
Shipping	Ice bag