

Anti-Clvs2 Rabbit Polyclonal Antibody For Wb, Elisa - Q5Syc1

Item Number: CM0023124

Introduction

Rabbit polyclonal antibody targeting human Clavesin-2 (CLVS2), validated for WB and ELISA. Ideal for studying endosomal/lysosomal morphology and neuronal protein trafficking. Recombinant immunogen.

[Learn More](#)

Parameter	Value
Product Name	CLVS2 Rabbit pAb
Remarks/Alias	RLBP1L2, C6orf212, C6orf213, bA160A10.4, CLVS2
Species	Human
Gene ID	134829
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 265-327 of human CLVS2 (NP_001010852.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q5SYC1
Protein Weight	38kDa
Shipping	Ice bag