

## Anti-Cmas Rabbit Polyclonal Antibody For Wb And Elisa - Q8Nfw8

Item Number: CM0015004

### Introduction

High-quality rabbit polyclonal antibody against human CMAS (CMP-N-acetylneuraminic acid synthase). Validated for WB and ELISA, cross-reacts with human, mouse, and rat. Ideal for sialic acid pathway research.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | CMAS Rabbit pAb                                                                                                  |
| Remarks/Alias    | CSS; CMAS                                                                                                        |
| Species          | Human                                                                                                            |
| GenelD (Human)   | 55907                                                                                                            |
| GenelD           | 55907                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-263 of human CMAS (NP_061156.1). |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, ELISA                                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | Q8NFW8                                                                                                           |
| Protein Weight   | 48kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |