

Anti-Cmip Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8ly22

Item Number: CM0019268

Introduction

High-quality rabbit polyclonal antibody against CMIP, validated for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for T-cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	CMIP Rabbit pAb
Remarks/Alias	TCMIP; CMIP
Species	Human
GeneID (Human)	80790
GeneID	80790
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 640-773 of human CMIP (NP_938204.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IY22
Protein Weight	86kDa
Shipping	Ice bag