


Anti-Cmklr1 Polyclonal Antibody For Wb And Elisa - Q99788

Item Number: CM0021350



Introduction

High-quality rabbit polyclonal antibody against human CMKLR1 (ChemR23), validated for WB and ELISA. Suitable for human, mouse, and rat samples. Ideal for metabolic inflammation and immune response research.

[Learn More](#)

Parameter	Value
Product Name	CMKLR1 Rabbit pAb
Remarks/Alias	DEZ; ERV1; RVER1; ChemR23; CHEMERINR; CMKLR1
Species	Human
GeneID (Human)	1240
GeneID	1240
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 250-350 of human CMKLR1 (NP_001135815.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99788
Protein Weight	42kDa
Shipping	Ice bag