


Anti-Cmklr1 Rabbit Polyclonal Antibody For Wb And Elisa - Q99788

Item Number: CM0014436



Introduction

High-quality rabbit polyclonal antibody against human CMKLR1 (Chemerin receptor) validated for WB and ELISA. Reacts with human and mouse. Ideal for immunology and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	CMKLR1 Rabbit pAb
Remarks/Alias	DEZ; ERV1; RVER1; ChemR23; CHEMERINR; CMKLR1
Species	Human
GeneID (Human)	1240
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 250-350 of human CMKLR1 (NP_001135815.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99788
Protein Weight	42kDa
Shipping	Ice bag