

Anti-Cmtm6 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Nx76

Item Number: CM0002386

Introduction

High-quality anti-CMTM6 rabbit monoclonal antibody for western blot (WB) and ELISA. Recognizes human and rat CMTM6, a master regulator of PD-L1 recycling and T-cell anergy. Ideal for immune checkpoint research and cancer immunology studies.

[Learn More](#)

Parameter	Value
Product Name	CMTM6 Rabbit mAb
Remarks/Alias	CKLFSF6; PRO2219; CMTM6
Species	Human
Gene ID (Human)	54918
Gene ID	54918
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 84-183 of human CMTM6 (NP_060271.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9NX76
Protein Weight	20kDa
Shipping	Ice bag