

Anti-Cnbp Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P62633

Item Number: CM0013267

Introduction

Anti-CNBP rabbit polyclonal antibody for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Suitable for studying transcriptional repression and RNA binding. 19 kDa protein.

[Learn More](#)

Parameter	Value
Product Name	CNBP Rabbit pAb
Remarks/Alias	DM2; ZNF9; CNBP1; PROMM; RNF163; ZCCHC22; CNBP
Species	Human
Gene ID	7555
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 70-170 of human CNBP (NP_001120668.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62633
Protein Weight	19kDa
Shipping	Ice bag