

Anti-Cndp1 Polyclonal Antibody For Wb, Elisa - Q96Kn2

Item Number: CM0030127



Introduction

CNDP1 Rabbit Polyclonal Antibody (CN1, CPGL2) validated for WB, ELISA. Targets human beta-Ala-His dipeptidase (serum carnosinase); cross-reacts with mouse, rat. For carnosine metabolism research.

[Learn More](#)

Parameter	Value
Product Name	CNDP1 Rabbit pAb
Remarks/Alias	CN1; CPGL2; HsT2308; CNDP1
Species	Human
Gene ID (Human)	84735
Gene ID	84735
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human CNDP1 (NP_116038.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96KN2
Protein Weight	57kDa
Shipping	Ice bag