

Anti-Cndp2 Rabbit Polyclonal Antibody For Wb And Elisa - Q96Kp4

Item Number: CM0013858



Introduction

Rabbit polyclonal antibody against CNDP2, validated for Western blot and ELISA in human, mouse, and rat. Cytosolic non-specific dipeptidase involved in glutathione metabolism and cell cycle regulation.

[Learn More](#)

Parameter	Value
Product Name	CNDP2 Rabbit pAb
Remarks/Alias	CN2; CPGL; PEPA; HsT2298; HEL-S-13; CNDP2
Species	Human
Gene ID (Human)	55748
Gene ID	55748
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human CNDP2 (NP_060705.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96KP4
Protein Weight	53kDa
Shipping	Ice bag