

Anti-Cnga2 Polyclonal Antibody For Wb, Elisa - Q16280

Item Number: CM0027294



Introduction

Rabbit polyclonal antibody against human CNGA2, validated for Western blot and ELISA. Reacts with human, mouse, rat. Ideal for olfactory signaling and cyclic nucleotide-gated ion channel research.

[Learn More](#)

Parameter	Value
Product Name	CNGA2 Rabbit pAb
Remarks/Alias	CNCA; CNG2; CNCA1; OCNC1; OCNCa; OCNCALPHA; CNGA2
Species	Human
GeneID (Human)	1260
GeneID	1260
Immunogen	Synthetic peptide; A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CNGA2 (NP_005131.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16280
Protein Weight	76kDa
Shipping	Ice bag