

Anti-Cnga3 Rabbit Polyclonal Antibody For Wb, Elisa - Q16281

Item Number: CM0027386



Introduction

Rabbit polyclonal antibody targeting human CNGA3, suitable for Western blot and ELISA. Recognizes cone photoreceptor cGMP-gated channel subunit alpha-3, key in phototransduction. Cross-reacts with Human, Mouse, Rat.

[Learn More](#)

Parameter	Value
Product Name	CNGA3 Rabbit pAb
Remarks/Alias	CNG3; ACHM2; CCNC1; CCNCa; CNCG3; CCNCalpha; CNGA3
Species	Human
GeneID (Human)	1261
GeneID	1261
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-165 of human CNGA3 (NP_001289.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16281
Protein Weight	79kDa
Shipping	Ice bag