

Anti-Cnga4 Rabbit Polyclonal Antibody For Wb, Elisa - Q8lv77

Item Number: CM0013013



Introduction

Rabbit polyclonal antibody against CNGA4, suitable for WB and ELISA applications, reacts with human and mouse samples. CNGA4 is a cyclic nucleotide-gated channel involved in olfactory signal transduction.

[Learn More](#)

Parameter	Value
Product Name	CNGA4 Rabbit pAb
Remarks/Alias	CNG4; CNG5; CNCA2; CNG-4; CNGB2; OCNC2; OCNCb; OCNCBETA; CNGA4
Species	Human
GeneID (Human)	1262
GeneID	1262
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 416-575 of human CNGA4 (NP_001032406.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8lv77
Protein Weight	66kDa
Shipping	Ice bag