

Anti-Cnmd Polyclonal Antibody For Wb, Elisa - 075829

Item Number: CM0023317

Introduction

Rabbit polyclonal antibody against human CNMD (CHM1/LECT1), validated for WB and ELISA. Regulates chondrocyte growth and angiogenesis; suitable for skeletal development and cardiovascular research.

[Learn More](#)

Parameter	Specification
Product Name	CNMD Rabbit pAb
Remarks/Alias	CHM1; CHM-I; LECT1; BRICD3; MYETS1; CNMD
Species	Human
GeneID (human)	11061
GeneID	11061
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 75-334 of human CNMD. (NP_008946.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O75829
Protein Weight	37kDa
Shipping	Ice bag