

Anti-Cnn3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q15417

Item Number: CM0024631

Introduction

Rabbit polyclonal antibody against human CNN3 (Calponin-3), suitable for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat CNN3 (~36 kDa). Ideal for smooth muscle contraction and actin-binding research.

[Learn More](#)

Parameter	Value
Product Name	CNN3 Rabbit pAb
Remarks/Alias	CNN3
Species	Human
Gene ID (Human)	1266
Gene ID	1266
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 220-329 of human CNN3 (NP_001830.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15417
Protein Weight	36kDa
Shipping	Ice bag