

Anti-Cnnm1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nru3

Item Number: CM0020112

Introduction

Rabbit polyclonal antibody against metal transporter CNNM1 (ACDP1, CLP-1). Validated for WB, IF/ICC, ELISA. Reacts with human, mouse, and rat. Immunogen: recombinant fusion protein covering amino acids 692-951 of human CNNM1.

[Learn More](#)

Parameter	Value
Product Name	CNNM1 Rabbit pAb
Remarks/Alias	ACDP1; CLP-1; CNNM1
Species	Human
GeneID (Human)	26507
GeneID	26507
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 692-951 of human CNNM1 (NP_065081.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NRU3
Protein Weight	104kDa
Shipping	Ice bag