


Anti-Cnnm4/Acdp4 Polyclonal Antibody For Wb, Elisa - Q6P4Q7

Item Number: CM0024569



Introduction

Rabbit polyclonal antibody targeting human CNNM4 (ACDP4), validated for WB and ELISA. Reacts with human, mouse, and rat. Ideal for research on metal transport, biomineralization, and retinal function.

[Learn More](#)

Parameter	Value
Product Name	CNNM4 Rabbit pAb
Remarks/Alias	ACDP4; CNNM4
Species	Human
Gene ID (Human)	26504
Gene ID	26504
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 700-775 of human CNNM4 (NP_064569.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6P4Q7
Protein Weight	87kDa
Shipping	Ice bag