

Anti-Cnot1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - A5Ykk6

Item Number: CM0030018



Introduction

High-quality rabbit polyclonal antibody against human CNOT1, validated for WB, IF/ICC, and ELISA. Recognizes human and mouse CNOT1. Ideal for investigating mRNA deadenylation and transcription regulation mechanisms.

[Learn More](#)

Parameter	Value
Product Name	CNOT1 Rabbit pAb
Remarks/Alias	NOT1; CDC39; HPE12; NOT1H; VIBOS; AD-005; CNOT1
Species	Human
Gene ID (Human)	23019
Gene ID	23019
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2137-2376 of human CNOT1 (NP_057368.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	A5YKK6
Protein Weight	267kDa
Shipping	Ice bag