


Anti-Cnot8 Polyclonal Antibody For Wb, Elisa - Q9Uff9

Item Number: CM0027527



Introduction

Anti-CNOT8 rabbit polyclonal antibody for WB and ELISA. Detects human and mouse CNOT8, a CCR4-NOT complex catalytic subunit involved in mRNA deadenylation and translational control. Suitable for molecular biology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	CNOT8 Rabbit pAb
Remarks/Alias	CAF1; POP2; CALIF; Caf1b; hCAF1; CNOT8
Species	Human
GeneID (Human)	9337
GeneID	9337
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-292 of human CNOT8 (NP_004770.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UFF9
Protein Weight	34kDa
Shipping	Ice bag