

Anti-Cnpase Polyclonal Antibody For Wb, Ihc-P, Elisa - P09543

Item Number: CM0019398



Introduction

Polyclonal antibody against human CNPase (CNP) validated for Western blot, IHC-P, and ELISA. Reacts with human, mouse, and rat samples. Ideal for neurobiology, myelin research, and antiviral studies.

[Learn More](#)

Parameter	Description
Product Name	CNPase Rabbit pAb
Remarks/Alias	CNP1; HLD20; CNPase
Species	Human
Gene ID	1267
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 152-421 of human CNPase (NP_149124.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09543
Protein Weight	48kDa
Shipping	Ice bag